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— Cautions:

- Do not bend, kink, or stretch the extension or adaptor, which may damage the component.

Cautions:

Do not bend, kink, or stretch the extension or adaptor, which may damage the component.

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- Place the pocket adaptor behind the neurostimulator. Coil the adaptor wire and the excess extension wire behind the adaptor, ensuring there are no sharp bends in any of the wires (Figure 10).

1. Place the pocket adaptor behind the neurostimulator. Coil the adaptor wire and the excess extension wire behind the adaptor, ensuring there are no sharp bends in any of the wires (Figure 10).

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
Figure 10. Ensure there are no sharp bends in any of the wires.

Sharp bend

There should be no sharp bends for the

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Тел.: (495) 580 73 77
Факс: (495) 580 73 78
www.medtronic.ru

adaptor and extension wires.

The adaptor and extension wires should not sharply bend.

YES

NO

Figure 10. Ensure there are no sharp bends in any of the wires.

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Caution: To use the nonsterile clinician programmer in a sterile field, place a sterile barrier between the patient and the programming head to prevent infection. Do not sterilize any part of the clinician programmer. Sterilization may damage the programmer.

Note: The neurostimulator should be in the pocket during system interrogation for integrity to ensure proper readings.

1. To ensure proper connection of each extension to the neurostimulator, use the clinician programmer to program the basic stimulation parameters, check the battery status, and check the electrode impedances to rule out a short or open circuit.
2. If the system integrity test results are not acceptable, refer to "Connecting the extension to the neurostimulator" on page 13.

Caution: To use the nonsterile clinical programmer in a sterile field, place a sterile barrier between the patient and the programming head to prevent infection. Do not sterilize any part of the clinical programmer. Sterilization may damage the programmer.

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1. To ensure proper connection of each extension to the neurostimulator, use the clinical programmer to program the basic stimulation parameters, check the battery

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status, and check the electrode impedances to rule out a short or open circuit.

If the system integrity test results are not acceptable, refer to "Connecting the extension to the neurostimulator" on page 13.

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