2569808

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
Instruction for Use

1. Product Description

Please read the instructions carefully before use!

The Angiosaver™ Inflation Device Kits contains an Inflation Device and a Y-Hemostasis Valve Set.
 Inflation Device is a kind of device which is used in PTCA operation for pressurizing on the balloon dilation catheter to make the balloon expand so as to realize the goal of expanding the blood vessel or releasing a stent inside the blood vessel.
 The Inflation Device is mainly consisted of a Syringe, a screwplunger, a reactor pad/Trigger, Shellupper and lower shell, a piston seat, a pressure gauge, a Tubing, pressure-resisting pipe, a male rotating adapter, cone locking joint and stopcock (optional).

Y-Hemostasis Valve Set is one of Disposable Interventional Accessories.

Y-Hemostasis Valve Set including a Y-Hemostasis Valve, a Guidewire Introducer and a Guidewire Torquer. Y-Hemostasis Valve set is mainly used for insertion of Interventional medical devices for Interventional diagnosis and treatment, hemostasis or liquid medicine injection. Guidewire Introducer can protect the Guidewire when it goes through. Guidewire Torquer is intended to fix the Guidewire and make it convenient for doctor to operate.

1. Y-Hemostasis Valve; 2. Guidewire Introducer; 3. Guidewire Torquer

2. Indications

Inflation Device Kits is indicated for use in percutaneous coronary intervention (PCI) treatment (such as percutaneous coronary angioplasty, PTCA, stent implantation operation) provides the pressure inside the balloon.

3. Contraindications

None.

4. Instructions

1. Before using, check the package and inflation devices to assure they are not damaged, please confirm if there is any damage caused by the package and the product.
2. Push forward the reactor pad/Trigger and pull backward the screwplunger at the same time. Release the reactor pad/Trigger with the dilation control medium (or other liquid) reaching the designated amount.
3. Put the Inflation Device straight up and push forward the reactor pad/Trigger, and pull backward the screwplunger at the same time to discharge the residual air within the Syringe. If necessary, release the reactor pad/Trigger and rotate the screwplunger appropriately to discharge the residual air within the Syringe.

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4. Connect the male rotating adapter/one locking joint of the Inflation Devices and the balloon catheter and keep the Pressure gauge dial plate facing directly towards the operator at the same time.
 5. Push forward the ejector pad/trigger and pull backward the screwplunger at the same time to evacuate the residual air inside the balloon catheter.
 6. Release the ejector pad/trigger, and rotate the screwplunger clockwise to pressurize the balloon; rotate the screwplunger anticlockwise to discharge the pressure of the balloon. When the pressure drops to a value lower than 10atm, push forward the ejector pad/trigger to discharge the pressure of the balloon. (Note: push forward the ejector pad/trigger so as to pull back or push forward the screwplunger freely, so the Inflation Devices will output negative pressure or positive pressure)
 7. Push forward the ejector pad/trigger, and pull backward the screwplunger at the same time. Release the ejector pad/trigger to form negative pressure inside the balloon.
 8. When the Inflation Devices is kept at a certain pressure, the pressure can be released quickly by pushing forward the ejector pad/trigger.
- Note: Over quick release of pressure might result in risks, please operate according to clinical conditions.

5. PRECAUTIONS:

- It is used only by professional doctors; read the instructions before using.
- This is a disposable device, repeated using is prohibited. The device is designed for single use. Please do not re-sterilize and/or reuse.
- Ethylene oxide sterilization is applied in this product; it belongs to a kind of sterile and non-pyrogenic product; so it is forbidden to be used if the package is damaged. Use prior to the USE-BY DATE.
- Improper fasten and connection Y-Hemostasis Valve Set with other devices may cause air bubble in the circulation system.
- Before using Inflation Devices, you must refer to the information related to the maximum bursting pressure specified by the manufacturer for balloon dilation catheters or other interventional devices. Pressure loss will indicate leakage of the system.
- Prevent Interventional Accessories from moisture, exposure, rain, high temperature, the weight and impact during transport.
- Disposing the device according to the local law after using.

6. MATERIALS TO BE USED IN COMBINATION WITH Y-HEMOSTASIS VALVE SET INCLUDED

- High Pressure Extension Tubing
- Manifold
- Guidewire
- Contrast medium

7. PREPARATION FOR USE

Examination before using: First check the package and don't use it if the package is damaged, then open the package to check the appearance and confirm that there is no defect such as loose of components, broken, and so on before using.


8. Storage:

It should be stored in the environment with the relative humidity lower than 80%, and the temperature lower than 40 C without any corrosive gas, and stored in a shady and cool, dry, well-ventilated and clean place.

9. Disclaimer of Warranty and Limitation of Remedy

There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the LEPU Medical product(s) described in this publication. Under no circumstances shall LEPU Medical be liable for any direct, incidental, or consequential damages resulting from the abnormal use of the product other than as expressly provided by specific law. No person has the authority to bind LEPU Medical to any representation or warranty except as specifically set forth herein.

10. Signs and legends:Symbol

Lepu Medical Technology (Beijing) Co., Ltd.

No. 37, Chaoyuan Road, Changping District, Beijing
102200, P.R. China
Tel: +86-10-80125566 Fax: +86-10-80120600

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